Table of Contents:

Symbols	4
Assembly of needle to syringe	6
Canada	8
English (en)	8
Français (fr)	10
Brazil	12
Português (pt)	12
Rest of the World	14
English (en)	14
Deutsche (de)	16
Français (fr)	18
Italiano (it)	20
Español (es)	22
Português (pt)	24
Русский (ru)	26

SYMBOL SYMBOL SYMBOLE SIMBOLO
 SIMBOLO SIMBOLOGIA СИМВОЛ:

Syringe fluid path sterilized using moist heat
 Sterilisation des Spritzenflüssigkeitspfades mittels feuchter Hitze
 Trajectoire du fluide dans la seringue stérilisée à la chaleur humide
 Percorso del fluido nella siringa sterilizzato a calore umido
 Canal del fluido de la jeringa esterilizado al vapor
 Esterilizado por vapor
 Путь протекающей по шприцу жидкости простерилизован с использованием влажного жара

Store between 2 and 25 °C
 Bei Temperaturen zwischen 2 °C und 25 °C lagern
 Conserver entre 2 et 25 °C
 Conservare ad una temperatura compresa tra 2°C e 25°C
 Almacenar entre 2 y 25°C
 Armazenar entre 2 e 25 °C
 Храните при температуре от 2°C до 25°C

Read the Instructions before using the product
 Lesen Sie die Anweisungen, bevor Sie das Produkt verwenden
 Lire les instructions avant d'utiliser le produit
 Leggere le istruzioni prima di utilizzare il prodotto
 Lea las instrucciones antes de usar el producto
 Leia as instruções de uso antes de usar o produto
 Перед использованием изделия прочтите все инструкции

Do not use if the package has been damaged
 Verwenden Sie das Produkt nicht, wenn die Verpackung beschädigt ist
 Ne pas utiliser si l'emballage a été endommagé
 Non utilizzare se la confezione è danneggiata
 No usar si el paquete se encuentra dañado
 Não use se a embalagem estiver aberta ou danificada
 Не используйте, если упаковка повреждена

Do not reuse
 Das Produkt darf nicht wiederverwendet werden
 Ne pas réutiliser
 Evitare il riutilizzo
 No reutilizar
 Produto de uso único - proibido reprocessar- não reutilizar
 Не используйте повторно

Keep away from sunlight
 Vor Sonneneinstrahlung schützen
 Tenir à l'écart de la lumière du soleil
 Non esporre alla luce solare
 Manterner protegido de la luz solar
 Não expora luz solar
 Не допускайте попадания солнечного света

Sterilized by irradiation
 Strahlensterilisation
 Stérilisé par irradiation
 Sterilizzato tramite irradiazione
 Esterilizado por irradiación
 Esterilizado por irradiação
 Стерилизован облучением

Lot number
 Chargennummer
 Numéro de lot
 Numero di lotto
 Número de LOTE
 Número do lote
 Номер партии

Date of Manufacture
 Herstellungsdatum
 Date de fabrication
 Data di produzione
 Fecha de fabricación
 Data de Fabricação
 Дата изготовления


Keep dry
 Trocken lagern
 Garder au sec
 Conservare in luogo secco
 Manterner seco
 Mantenha seco
 Храните в сухом месте

Expiration date
 Ablaufdatum
 Date d'expiration
 Data di scadenza
 Fecha de caducidad
 Data de validade
 Срок годности

Manufacturer
 Hersteller
 Fabricant
 Produttore
 Fabricante
 Fabricante
 Производитель

ASSEMBLY OF NEEDLE TO SYRINGE:

1. Carefully unscrew the syringe end cap.
2. Firmly grip the wider part of the needle guard.
3. Press and rotate the needle guard by 90° (one quarter turn). This quarter turn is necessary to lock the needle into the syringe. Remove the needle guard.
4. Using the thumb and forefinger, tightly secure the glass syringe cylinder and the Luer-lock adapter.
5. Grip the needle guard with the opposite hand. To facilitate proper assembly, push and rotate firmly.

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
DESCRIPTION

Revanesse® Outline™ is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1.2ml syringes of Revanesse® Outline™ along with two sterilized 27G needles.

APPLICATION RANGE / INDICATIONS

Revanesse® Outline™ is a cross linked hyaluronic acid gel that is indicated for the treatment of medium to deep facial rhytides and lipotrophy by injection into the middle part of the dermis layer. Implantation site is dependent on depth and location of injection, and averages 6-12 months.

ANTICIPATED SIDE EFFECTS

Physicians must inform patients that with every injection of Revanesse® Outline™ there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.

• Glabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis.

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CONTRAINDICATIONS

- Do not inject Revanesse® Outline™ into eye contours (into the eye circle or eyelids).
- Pregnant women or lactating women should not be treated with Revanesse® Outline™.
- Revanesse® Outline™ is only intended for intradermal use and must not be injected into blood vessels. This may occlude and could cause an embolism.
- Patients who develop hypertrophic scarring should not be treated with Revanesse® Outline™.
- Patients that have a known history of streptococcal disease.
- Never use Revanesse® Outline™ in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.
- People under the age of 18 should not be treated with Revanesse® Outline™.
- Patients with acne and / or other inflammatory diseases of the skin should not be treated with Revanesse® Outline™.
- Patients with unattainable expectations.
- Patients with multiple severe allergies.
- Patients with acute or chronic skin disease in or near the injection sites.
- Coagulation defects or under anti-coagulation therapy.
- Patients with sensitivity to hyaluronic acid.

It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (ie: corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revanesse® family of products and / or to Prollemium Medical Technologies Inc. directly.

ADMINISTRATION & DOSAGE

- Revanesse® Outline™ should only be injected by or under the direct supervision of qualified physicians who have been trained on the proper injection technique for filling facial wrinkles.
- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.
- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.
- Inject the product slowly and apply the least amount of pressure necessary.
- Revanesse® Outline™ and needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.
- Keep the product at room temperature for 30 minutes prior to injection.
- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be massaged until the skin returns to its normal color.
- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.

PRECAUTIONS

- Revanesse® Outline™ should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.

- Revanesse® Outline™ should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revanesse® Outline™ never comes into contact with this substance or medical instrumentation that has come into contact with this substance.

- Revanesse® Outline™ should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.
- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.
- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.
- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.
- If you are taking aspirin, non-steroidal anti-inflammatory medications, St. John's Wort, High doses of Vitamin E supplements or any similar medications prior to treatment, be aware that these may increase bruising and bleeding at the injection site.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

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- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.
- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- Physicians should take extra precautions when treating patients with auto immune disorder, HIV, AIDS or who are under immuno-therapy as these patients are at higher risk of adverse events.
- **NOTE:** The correct injection technique is crucial to treatment success & patient satisfaction. Revanesse® Outline™ should only be injected by a practitioner qualified according to local laws and standards.
- The graduation on the syringe is not precise and should be used as a guide only. The amount of material to be injected is best determined by visual and tactile assessment by the user.

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FABRICANTE

Prolenium Medical Technologies, Inc.
138 Industrial Parkway North, Aurora, ON
L4G 4C3, Canada

T: 1.866.353.3015 | F: 1.866.876.6716
T: 905.508.1469 (local) | Fax: 905.508.6716 (local)
E: info@prolenium.com www.prolenium.com

COMPOSITION

Cross-linked hyaluronic acid (High Viscosity) 25mg/ml
 In phosphate buffered saline
 [Cross-linked with Butanediol-diglycidylether (BDDE)]

DESCRIPTION

Revasse® Outline™ is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1.2ml syringes of Revasse® Outline™ along with two sterilized needles.

APPLICATION RANGE / INDICATIONS

Medical Indications: The products are space-occupying tissue reconstructive materials composed of a hyaluronic acid gel that is indicated for restoration of volume lost from lipatrophy / lipodystrophy, and / or correction of contour deficiencies and anatomic deformities of either pathologic origin or after trauma, in soft tissue.

Intended patients are those desiring the correction of contour deficiencies and deformities of soft tissue, such as HIV-associated lipatrophy and lipodystrophy.

Cosmetic Indication: Revasse® Outline™ is indicated for the treatment of facial rhytides, volume restoration, lip augmentation, skin hydration and contouring of depressions by injection into tissue.

ANTICIPATED SIDE EFFECTS

Physicians must inform patients that with every injection of Revasse® Outline™ there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.

- Nodules or induration are also possible at the injection site.

- Poor product performance due to improper injection technique.

- Glabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis.

Reactions thought to be of hypersensitivity in nature have been reported in less than one in every 1500 treatments. These have consisted of prolonged erythema, swelling and induration at the implant site.

These reactions have started either shortly after injection or after a delay of 2 – 4 weeks and have been described as mild or moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

CONTRAINDICATIONS

- Do not inject Revasse® Outline™ into eye contours (into the eye circle or eyelids).

- Pregnant women or lactating women should not be treated with Revasse® Outline™.

- Revasse® Outline™ is only intended for intradermal use and must not be injected into blood vessels. This may occlude and could cause an embolism.

- Patients who develop hypertrophic scarring should not be treated with Revasse® Outline™.

- Patients that have a known history of streptococcal disease.

- Never use Revasse® Outline™ in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.

- People under the age of 18 should not be treated with Revasse® Outline™.

- Patients with acne and / or other inflammatory diseases of the skin should not be treated with Revasse® Outline™.

- Patients with unattainable expectations.

- Patients with multiple severe allergies.

- Patients with acute or chronic skin disease in or near the injection sites.

- Coagulation defects or under anti-coagulation therapy.

- Patients with sensitivity to hyaluronic acid.

It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (ie: corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revasse® family of products and / or to Prollemium Medical Technologies Inc. directly.

ADMINISTRATION & DOSAGE

- Revasse® Outline™ should only be injected by or under the direct supervision of qualified physicians who have been trained on the proper injection technique for filling facial wrinkles.

- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.

- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.

- Inject the product slowly and apply the least amount of pressure necessary.

- Revasse® Outline™ needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.

- Keep the product at room temperature for 30 minutes prior to injection.

- The depth of injection is at the discretion of the practitioner/physician and is dependent on the severity of the wrinkle/depression and the level of correction desired.

- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be massaged until the skin returns to its normal color.

- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.

PRECAUTIONS

- Revasse® Outline™ should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.

- Revasse® Outline™ should not be injected into an area where there is a permanent filler or implant.

- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revasse® Outline™ never comes into contact with this substance or medical instrumentation that has come into contact with this substance.

- Revasse® Outline™ should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.

- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.

- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.

- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.

- If you are taking aspirin, non-steroidal anti-inflammatory medications, St. John's Wort, High doses of Vitamin E supplements or any similar medications prior to treatment, be aware that these may increase bruising and bleeding at the injection site.

- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

WARNINGS

Confirm that the seal on the box has not been broken and sterility has not been compromised. Confirm that the product has not expired. Product is for single use only; do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.

- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.

- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

Physicians should take extra precautions when treating patients with auto immune disorder, HIV, AIDS or who are under immuno-therapy as these patients are at higher risk of adverse events.

NOTE: The correct injection technique is crucial to treatment success & patient satisfaction. Revasse® Outline™ should only be injected by a practitioner qualified according to local laws and standards.

The graduation on the syringe is not precise and should be used as a guide only. The amount of material to be injected is best determined by visual and tactile assessment by the user.

SHELF LIFE & STORAGE

Expiry is indicated on each individual package. Store between 2°-25° C, and protect from direct sun light and freezing.

THE MANUFACTURER

Prollemium Medical Technologies, Inc.
 138 Industrial Parkway North, Aurora, ON
 L4G 4C3, Canada

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Prolium Medical Technologies, Inc.