The Revanesse® Family of Products

ReDexis® Ultra^{W710}
Cross Linked Hyaluronic Acid
with Lidocaine

Dextranomers +

SYRINGE

Syringe fluid path sterilized using moist heat
Fluidweg in der Spritze mit feuchter Hitze sterilisiert
Trajet de fluide dans la seringue stérilisée à la chaleur humide
Canale del fluido nella siringa sterilizzato a calore umido
Trayectoria de fluido en la jeringa esterilizada usando calor húmedo
Passagem do fluido na seringa esterilizada utilizando calor úmido
Содержимое шприца было простерилизовано влажным жаром

Do not use if the package has been damaged
Nicht verwenden, falls die Verpackung beschädigt ist
Ne pas utiliser si l'emballage a été endommagé
Non utilizzare se la confezione è stata danneggiata
No usar si el paquete está dañado
Não use se o pacote estiver danificado
Не использовать, если упаковка была повреждена

NEEDLE

Revitalize. Rejuvenate. ReNew®

COMPOSITION

Cross-linked hyaluronic acid	17mg/ml
DEAE Sephadex A25	50mg/ml
Lidocaine	0.3% w/w
In phosphate buffered saline	
[Cross linked with Butandiol-diglycidylether (BDDE)]	

DESCRIPTION

ReDexis® Ultra + is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1.2ml syringes of ReDexis® Ultra + along with two sterilized needles.

APPLICATION RANGE / INDICATIONS

ReDexis® Ultra™ + is a cross linked hyaluronic acid gel that is indicated for the treatment of medium to deep facial rhytides and lipotrophy by injection into the middle part of the dermis layer. Implantation life is dependent on depth and location of injection, and averages 6-12 months.

ANTICIPATED SIDE EFFECTS

Physicians must inform patients that with every injection of ReDexis® Ultra + there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.
- Glabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis.

Reactions thought to be of hypersensitivity in nature have been reported in less than one in every 1500 treatments. These have consisted of prolonged erythema, swelling and induration at the implant site.

These reactions have started either shortly after injection or after a delay of 2 – 4 weeks and have been described as mild or moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

CONTRAINDICATIONS

- Contains lidocaine and is contraindicated for patients with a history of allergies to such material.
- Do not inject ReDexis® Ultra + into eye contours (into the eye circle or eyelids).
- Pregnant women, or women during lactation should not be treated with ReDexis® Ultra +.
- ReDexis® Ultra + is only intended for intradermal use and must not be injected into blood vessels. This may occlude and could cause an embolism.
- Patients who develop hypertrophic scarring should not be treated with ReDexis® Ultra +.
- Contains trace amounts of gram positive bacterial proteins, and is contraindicated for patients with a history of allergies to such material.
- Never use ReDexis® Ultra + in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.
- People under the age of 18 should not be treated with ReDexis® Ultra +.
- Patients with acne and / or other inflammatory diseases of the skin should not be treated with ReDexis® Ultra +.
- Patients with unattainable expectations.
- Patients with auto-immune disorders or under immunotherapy.
- Patients with multiple severe allergies.
- Patients with acute or chronic skin disease in or near the injection sites.
- Coagulation defects or under anti-coagulation therapy.
- Patients with sensitivity to hyaluronic acid.

It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (ie: corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revanesse® family of products and / or to Prolenium Medical Technologies Inc. directly.

ADMINISTRATION & DOSAGE

- ReDexis® Ultra + should only be injected by or under the direct supervision of qualified physicians who have been trained on the proper injection technique for filling facial wrinkles.
- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.
- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.
- Inject the product slowly and apply the least amount of pressure necessary.
- ReDexis® Ultra + and needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.
- Keep the product at room temperature for 30 minutes prior to injection.
- The amount of ReDexis® Ultra + that is injected will depend on the severity of the depression or the level of correction desired.
- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be

massaged until the skin returns to its normal color.

- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.
- ReDexis® Ultra + must be injected into the deep dermis.

PRECAUTIONS

- ReDexis® Ultra + should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.
- ReDexis® Ultra + should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that ReDexis® Ultra + never comes into contact with this substance or medical instrumentation that has come into contact with this substance.
- ReDexis® Ultra + should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.
- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.
- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.
- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.
- If you are using aspirin, non-steroidal anti-inflammatory medications, St. John's Wort or High doses of Vitamin E supplements prior to treatment or any similar medications be aware that these may increase bruising and bleeding at the injection site.
- Based on a toxicological risk assessment for lidocaine, patients should be limited to 20 mL per 60 kg (130 lbs) body mass per annum. The safety of injecting greater amounts has not been established.
- The safety for use in patients under 18 years or over 65 years has not been established.
- Patients who are visibly ill, with bacterial or viral infections, influenza, or active fever should not be treated until a resolution of their symptoms.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

WARNINGS

Confirm that the seal on the box has not been broken and sterility has not been compromised. Confirm that the product has not expired. Product is for single use only; do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.

- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.
- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

The Revanesse family of products should not be used in areas that have high vascularity. Use in these areas such as the glabella and nose region has resulted in cases of vascular embolization and symptoms consistent with ocular vessel occlusion (i.e.: blindness).

SHELF LIFE & STORAGE

Expiry is indicated on each individual package. Store between 2°-25° C, and protect from direct sun light and freezing.

NOTE: The correct injection technique is crucial to treatment success & patient satisfaction. ReDexis® Ultra + should only be injected by a practitioner qualified according to local laws and standards.

The graduation on the syringe is not precise and should be used as a guide only. The amount of material to be injected is best determined by visual and tactile assessment by the user.

MANUFACTURER

Prollenium Medical Technologies, Inc.
138 Industrial Parkway North, Aurora, ON
L4G 4C3, Canada

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