---

REVANESSE®  
KISS™+

---

Instructions for Use  
Mode d'emploi

## Table of Contents:

Symbols .....	4
Assembly of needle to syringe.....	5
Canada .....	6
English (en).....	6
Français (fr) .....	8

## SYMBOL SYMBOLE

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Date of Manufacture  
Date de fabrication



Keep dry  
Garder au sec

Expiration date  
Date d'expiration

Manufacturer  
Fabricant

### ASSEMBLY OF NEEDLE TO SYRINGE:

1. Carefully unscrew the syringe end cap.
2. Firmly grip the wider part of the needle guard.
3. Press and rotate the needle guard by 90° (one quarter turn). This quarter turn is necessary to lock the needle into the syringe. Remove the needle guard.
4. Using the thumb and forefinger, tightly secure the glass syringe cylinder and the Luer-lock adapter.
5. Grip the needle guard with the opposite hand. To facilitate proper assembly, push and rotate firmly.

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**COMPOSITION**

Cross-linked hyaluronic acid .....	25mg/ml
Lidocaine .....	0.3% w/w
In phosphate buffered saline	
[Cross linked with Butandiol-diglycidylether (BDDE)]	

**DESCRIPTION**

Revanesse® Kiss™+ is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1ml syringes of Revanesse® Kiss™+ along with two sterilized needles.

**APPLICATION RANGE / INDICATIONS**

Revanesse® Kiss™+ is a cross linked hyaluronic acid gel that is indicated for submucosal implantation for lip augmentation.

**ANTICIPATED SIDE EFFECTS**

Physicians must inform patients that with every injection of Revanesse® Kiss™+ there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.
- Glabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis. Reactions thought to be of hypersensitivity in nature have been reported in less than one in every 1500 treatments. These have consisted of prolonged erythema, swelling and induration at the implant site.

These reactions have started either shortly after injection or after a delay of 2 – 4 weeks and have been described as mild or moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

**CONTRAINDICATIONS**

- Revanesse® Kiss™+ is only intended submucosal injection into the lips and must not be injected into blood vessels. Implantation of Revanesse® Kiss™+ into dermal vessels may cause vascular occlusion, infarction, or embolic phenomena.
- Do not inject Revanesse® Kiss™+ into eye contours (into the eye circle or eyelids).
- Pregnant women, or women during lactation should not be treated with Revanesse® Kiss™+.
- Patients who develop hypertrophic scarring or keloid formation should not be treated with Revanesse® Kiss™+.
- Patients with evidence of scars at the intended treatment sites should not be treated with Revanesse® Kiss™+.
- Patients with known susceptibility to hypigmentation.
- Contains trace amounts of gram positive bacterial proteins, and is contraindicated for patients with a history of allergies to such material.
- Patients with allergies to hyaluronic acid products, or Streptococcal proteins should not be treated with Revanesse® Kiss™+.
- Revanesse® Kiss™+ contains lidocaine, and is contraindicated for patients with a history of allergies or sensitivities to such material and should not be used in patients with previous hypersensitivity to local anesthetics of the amide type, such as lidocaine.

- Never use Revanesse® Kiss™+ in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments, with Over-the-counter (OTC) wrinkle products or prescription wrinkle treatments within 4 weeks (28 days) prior to treatment.
- People under the age of 18 should not be treated with Revanesse® Kiss™+.
- Patients with acne and / or other inflammatory diseases of the skin should not be treated with Revanesse® Kiss™+.
- Patients with unattainable expectations.
- Patients with auto-immune disorders or under immunotherapy.
- Patients who have plans to undergo desensitization therapy should not be treated with Revanesse® Kiss™+.
- Patients with multiple severe allergies, or with allergy history including anaphylaxis, multiple severe allergies, atopy, should not be treated with Revanesse® Kiss™+.
- Patients with acute or chronic skin disease in or near the injection sites, or with any infection or unhealed wound of the face.
- Patients who are under concomitant anticoagulant therapy, antiplatelet therapy, or history of bleeding disorders, coagulation defects or connective tissue disorders should not use this product.

It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (ie: corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor or the Revanesse® family of products and / or ProLium Medical Technologies Inc. directly.

**ADMINISTRATION & DOSAGE**

- Revanesse® Kiss™+ should only be injected by or under the direct supervision of qualified physicians who have been trained in the proper injection technique for filling facial wrinkles.
- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.
- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.

- Inject the product slowly and apply the least amount of pressure necessary.
- Revanesse® Kiss™+ and needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.
- Keep the product at room temperature for 30 minutes prior to injection.
- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be massaged until the skin returns to its normal color.
- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.

**PRECAUTIONS**

- Revanesse® Kiss™+ should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.
- Revanesse® Kiss™+ should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revanesse® Kiss™+ never comes into contact with this substance or medical instrumentation that has come into contact with this substance.
- Revanesse® Kiss™+ should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.
- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.
- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.
- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.
- If you are using aspirin, non-steroidal anti-inflammatory medications, St. John's Wort of High doses of Vitamin E supplements prior to treatment or any similar medications be aware that these may increase bruising and bleeding at the injection site.
- Based on a toxicological risk assessment for lidocaine, patients should be limited to 20 mL per 60 kg (130 lbs) body mass per annum. The safety of injecting greater amounts has not been established.
- The safety for use in patients under 18 years or over 65 years has not been established.
- Patients who are visibly ill, with bacterial or viral infections, influenza, or active fever should not be treated until a resolution of their symptoms.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

**WARNINGS**

Confirm that the seal on the box has not been broken and sterility has not been compromised. Confirm that the product has not expired. Product is for single use only; do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.

- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include stroke, skin necrosis, and damage to underlying facial structures.
- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

The Revanesse® family of products should not be used in areas that have high vascularity. Use in these areas such as the glabella and nose region has resulted

**SHELF LIFE & STORAGE**

Expiry is indicated on each individual package. Store between 2°-25° C, and protect from direct sun light and freezing.

**NOTE:** The correct injection technique is crucial to treatment success & patient satisfaction. Revanesse® Kiss™+ should only be injected by a practitioner qualified according to local laws and standards.

The graduation on the syringe is not precise and should be used as a guide only. The amount of material to be injected is best determined by visual and tactile assessment by the user.

**MANUFACTURER**

ProLium Medical Technologies, Inc.  
138 Industrial Parkway North, Aurora, ON  
L4G 4C3, Canada

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Prolenium Medical Technologies, Inc.